

BONE WAX

BONE WAX

Advantages

- Soft and easy to form
- Easy to use
- Double packaged
- Gamma sterilized



Our Operation Cares About Your Operation...

CE 2764

BONE WAX

BONE WAX

Bonewax is a sterile haemostatic material to control of bleeding in bone injuries.

Hemostatic effect will be produce mechanically by occlusion of intaosseous vessels during surgical procedures. Bonewax is produced from PH Eur beeswax.It is recomended to knead or mould the Bonewax when it is "hand warm" to shape it before using.

Indications

- In orthopedic sugery & traumatology
- In thoracic surgery – sternum and costa
- In dental & maxillofacial surgery
- In neurosurgery – trepanation
- In general and plastic surgery



Type	Size	Pcs/Box	Reference
Bonewax	2,5 gr	12	BW12
Bonewax	2,5gr	24	BW24



 Yücel Medikal ve Tekstil Ürünleri San. Tic. Ltd. Şti.
Selahaddin Eyyubi Mah.1612 Sok. No:19
34850 Esenyurt / İstanbul / TURKEY
+90 212 672 86 10
+90 212 672 86 09
+90 544 672 86 09
info@yucelmedikal.com
www.yucelmedikal.com

CE 2764