

CLINISPONGE®

Absorbable Haemostatic Gelatin Sponge

Clinisponge® is an implantable, absorbable gelatin sponge with haemostatic effect and used in surgical procedures to achieve fast and effective hemostasis.


The sponge has a uniform porosity and reacts naturally to the coagulation process.

Properties & Advantages

- Biocompatible and completely resorbable
- High absorption capacity - up to 35 times of its own weight
- Resorbed approx. in 3-4 weeks
- Used as a carrier for drugs
- Used dry or soaked with sterile sodium chloride solution
- Well adhering ability to the bleeding site
- Variety of sizes
- Can be cut to desired size




CE 2460



Our Operation Cares About Your Operation...

CLINISPONGE®

Absorbable Haemostatic Gelatin Sponge



Indications

- General surgery
- Orthopaedic surgery
- Thoracic surgery
- Neurosurgery
- Gynecology
- Urology
- Maxillofacial surgery
- Orthopaedic surgery
- Thoracic surgery

Type	Size	Pcs/Box	Reference
Clinisponge Standart	80 x 50 x 10 mm	10	CL5810
Clinisponge Standart	70 x 50 x 10 mm	10	CL5710
Clinisponge Special	80 x 50 x 1 mm	10	CL5811
Clinisponge Special	70 x 50 x 1 mm	10	CL5711
Clinisponge Film	200 x 70 x 0,5 mm	10	CL2075
Clinisponge Dial	30 x 30 x 10 mm	10	CL3010
Clinisponge Tampon	80 x ø30 mm	5	CL8030
Clinisponge Dental	10 x 10 x 10 mm	50	CL1010
Clinisponge Dental	14 x 7 x 7 mm	50	CL1477
Clinisponge Powder	1 gr	5	CL1000

CE 2460



Yücel Medikal ve Tekstil Ürünleri San. Tic. Ltd. Şti.
Selahaddin Eyyubi Mah.1612 Sok. No:19
34850 Esenyurt / İstanbul / TURKEY
+90 212 672 86 10
+90 212 672 86 09
+90 544 672 86 09
info@yucelmedikal.com



www.yucelmedikal.com

