

Surgical Local Hemostats



CLINISPONGE®

Absorbable Haemostatic Gelatin Sponge



LYOCOLL®

Absorbable Collagen Haemostatic Felt & Powder



BONE WAX

Cliniwax 2,5 g

CLINISPONGE®

Absorbable Haemostatic Gelatin Sponge

Clinisponge® is an implantable, absorbable gelatin sponge with haemostatic effect and used in surgical procedures to achieve fast and effective hemostasis. The sponge has a uniform porosity and reacts naturally to coagulation process.

Properties & Advantages

- Biocompatible and completely resorbable
- High absorption capacity - up to 50 times of its own weight
- Resorbed approx. in 3-4 weeks
- Used as a carrier for drugs
- Used dry or soaked with physiological solution
- Well adhering ability to the bleeding site
- Variety of sizes
- Can be cut to desired size

Indications

- General surgery
- Orthopaedic surgery
- Thoracic surgery
- Neurosurgery
- Gynecology
- Urology
- Maxillofacial surgery
- Orthopaedic surgery
- Thoracic surgery



Type	Size	Pcs/Box	Reference
Clinisponge Standart	80 x 50 x 10 mm	10	CL5810
Clinisponge Standart	70 x 50 x 10 mm	10	CL5710
Clinisponge Special	80 x 50 x 1 mm	10	CL5811
Clinisponge Special	70 x 50 x 1 mm	10	CL5711
Clinisponge Film	200 x 70 x 1,5 mm	10	CL2075
Clinisponge Dial	30 x 30 x 10 mm	10	CL3010
Clinisponge Tampon	80 x ø30 mm	5	CL8030
Clinisponge Dental	10 x 10 x 10 mm	50	CL1010
Clinisponge Dental	14 x 7 x 7 mm	50	CL1477
Clinisponge Powder	1 gr	5	CL1000

Custom sizes may be produced according to your orders.

LYOCOLL®

Absorbable Collagen Haemostatic Felt & Powder


Lyocoll is a nativ collagen haemostatic material in fibrillar structure in bovine origin and practically used in all types of surgery to control of capillary, venous bleeding and oozing.

Lyocoll has an excellent haemostatic property and its surface structure also provides an additional strength to the coagulation procedure.

Lyocoll Collagen products are more effectively for local haemostasis to compareit with all haemostyptics made from other basic materials.

Properties & Advantages

- Excellent haemostatic effect
- Provides additional strength to coagulation
- Indicated for all surgical procedures
- Suitable for use on flat surfaces and to wrap vessels and anastomosis sites
- High tissue compatibility
- Excellent adherence to bleeding surface
- Absorbed within 3 weeks



Type	Size	Pcs/Box	Reference
Lyocoll Felt	3 x 5 cm	10	LY35
Lyocoll Felt	5 x 8 cm	10	LY58
Lyocoll Felt	10 x 12 cm	5	LY1012
Lyocoll Powder	0,5 gr	5	LYP05
Lyocoll Powder	1,0 gr	5	LYP10
Lyocoll Powder	5,0 gr	1	LYP50

Custom sizes may be produced according to your orders.

BONE WAX

Cliniwax 2,5 g

Cliniwax is a sterile haemostatic material to control of bleeding in bone injuries. Hemostatic effect will be produce mechanically by occlusion of intaosseous vessels during surgical procedures. Cliniwax is produced from PH Eur beeswax. It is recomended to knead or mould the Cliniwax when it is "hand warm" to shape it before using.

Indications

- In orthopedic sugery & traumatology
- In thoracic surgery – sternum and costa
- In dental & maxillofacial surgery
- In neurosurgery – trepanation
- In general and plastic surgery

Advantages

- Soft and easy to form
- Easy to use
- Double packaged
- Gamma sterilized



Type	Size	Pcs/Box	Reference
Bonewax Cliniwax	2,5 gr	12	BW12
Bonewax Cliniwax	2,5gr	24	BW24

Custom sizes may be produced according to your orders.

CE 1014



YÜCEL MEDİKAL Ltd. Şti.
Örnek Mah. 1612 Sok. Evren 2 Sanayi Sitesi
24. Blok No:19 34850 Esenyurt
İstanbul / TURKEY
+90 212 672 86 10
+90 212 672 86 09
+90 544 672 86 09
info@yucelmedikal.com



www.yucelmedikal.com