

# BONE WAX

## Cliniwax 2,5 g

## BONE WAX

### Advantages

- Soft and easy to form
- Easy to use
- Double packaged
- Gamma sterilized



# BONE WAX

## Cliniwax 2,5 g

## BONE WAX

**Cliniwax is a sterile haemostatic material to control of bleeding in bone injuries.**

Hemostatic effect will be produce mechanically by occlusion of intaosseous vessels during surgical procedures. Cliniwax is produced from PH Eur beeswax.It is recomended to knead or mould the Cliniwax when it is "hand warm" to shape it before using.

### Indications

- In orthopedic sugery & traumatology
- In thoracic surgery – sternum and costa
- In dental & maxillofacial surgery
- In neurosurgery – trepanation
- In general and plastic surgery



Type	Size	Pcs/Box	Reference
Bonewax Cliniwax	2,5 gr	12	BW12
Bonewax Cliniwax	2,5gr	24	BW24